
OBJECTIVE

Measuring the period of oscillation of a string pendulum with bobs of various masses.

EXPERIMENT PROCEDURE

- Measure the period of oscillation T of a string pendulum as a function of the length of the pendulum L .
- Measure the period of oscillation T of a string pendulum as a function of the mass of the pendulum bob m .
- Determine the acceleration due to gravity g .

SUMMARY

The period of oscillation T for a string pendulum is dependent on the length of the pendulum L , but does not depend on the mass of the bob m . This is to be verified by a series of measurements in which the period of oscillation of such a pendulum is measured by means of a photoelectric sensor connected to a digital counter.

REQUIRED APPARATUS

Quantity	Description	Number
1	Set of 4 Pendulum Bobs	1003230
1	Cord for Experiments	1001055
1	Tripod Stand 185 mm	1002836
1	Stainless Steel Rod 1500 mm	1002937
1	Stainless Steel Rod 100 mm	1002932
1	Clamp with Hook	1002828
2	Universal Clamp	1002830
1	Photo Gate	1000563
1	Digital Counter (230 V, 50/60 Hz)	1001033 or
	Digital Counter (115 V, 50/60 Hz)	1001032
1	Pocket Measuring Tape, 2 m	1002603
1	Electronic Scale 200 g	1003433

1
BASIC PRINCIPLES

A string pendulum with a bob of mass m and a length L will exhibit simple harmonic oscillation about its rest point as long as the angle of deflection is not too great. The period T , i.e. the time it takes for the pendulum to swing from one end of its motion to the other end and back, is dependent solely on the length of the pendulum L and not on the mass m .

If the pendulum is deflected from its rest position by an angle φ , the restoring force is as follows:

$$(1a) \quad F_1 = -m \cdot g \cdot \sin \varphi.$$

For small angles φ , this closely approximates to the following:

$$(1b) \quad F_1 = -m \cdot g \cdot \varphi$$

The moment of inertia of the accelerated mass is given by

$$(2) \quad F_2 = m \cdot L \cdot \ddot{\varphi}$$

Both these forces are equal, thus the result is equivalent for the equation of motion for simple harmonic oscillation:

$$(3) \quad \ddot{\varphi} + \frac{g}{L} \cdot \varphi = 0$$



For the period of oscillation T the following applies:

$$(4) \quad T = 2\pi \cdot \sqrt{\frac{L}{g}}$$

In this experiment the period of oscillation will be measured for various lengths of pendulum and masses of bob with the help of a photoelectric sensor connected to a digital counter. The digital counter's internal programming is such that it halts the time measurement after each complete swing of the pendulum.

EVALUATION

The measurements are plotted on a graph of T against L and another one of T against m . These graphs will show that the period of oscillation depends on the pendulum's length and not on the mass of the bob, as expected.


 Fig. 1: Period of oscillation T as a function of the pendulum length L

 Fig. 2: Period of oscillation T as a function of the pendulum mass m