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Safety and Warning Instructions for Kinesiology Tapes

General Guidelines

4. **Consult a Professional:** If you are uncertain or have health concerns, consult a doctor, physiotherapist, or qualified professional before using kinesiology tapes.
5. **Use Only as Intended:** The tape is designed for external use only. Do not apply to mucous membranes or open wounds.
6. **Skin Preparation:** Ensure the skin is clean, dry, and free of oils, creams, or dirt to ensure optimal adhesion and effectiveness.

Warnings

1. **Skin Reactions**
 - If irritation, redness, or itching occurs, remove the tape immediately.
 - Individuals with sensitive skin or adhesive allergies should conduct a patch test before use.
2. **Do Not Apply to Open Wounds:** Avoid using the tape on open wounds, irritated skin, burns, or skin infections.
3. **Circulation Issues:** Individuals with circulatory disorders or lymphatic problems should only use the tape after consulting a doctor.

4. **Contraindications:**

- Do not use in cases of thrombosis, deep vein thrombosis (DVT), or known blood clotting disorders.
- Avoid use with heart failure, skin cancer, or other severe medical conditions without prior medical approval.

5. **For Children and Pregnant Individuals:**

- Use on children or during pregnancy only after consulting a healthcare professional.
- Exercise caution, especially when applying tape to the abdomen or lower back during pregnancy.

6. **Do Not Apply Too Tightly:** Overly tight application can impair circulation. Light tension is usually sufficient for the desired effect.

7. **Prolonged Use:** Do not leave the tape on the skin longer than the recommended duration (usually 3-5 days). Allow the skin to recover before reapplying.

8. **Temperature and Moisture Exposure:** Avoid extreme heat, cold, or excessive sweating, as these can impact the tape's adhesion and function.

Additional Safety Tips

3. **Tape Removal:**

- Remove the tape gently to avoid skin irritation. Moistening the tape with oil or warm water can help.
- Avoid pulling the tape off abruptly to prevent pain or skin damage.

4. **Single-Use Only:** Kinesiology tapes are single-use products. Reusing them can compromise adhesion and hygiene.

Consult a Doctor if Symptoms Occur

If you experience pain, swelling, or other unusual symptoms, remove the tape immediately and seek medical advice.

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