1026563 Diodenlaser K, 532 nm, grün

DE

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09/25 HJB

LASER SICHERHEIT – Klasse 2 (IEC/EN 60825-1), 532 nm, max. 1 mW

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LASER SAFETY – Class 2 (IEC/EN 60825-1), 532 nm, max. 1 mW

General Safety

1. The laser emits visible radiation with a wavelength of 532 nm at a maximum output power of less than 1 mW, thus complies to class 2 regulations as specified in DIN EN 60825-1 "Safety of lasers", i.e. the human eye can be protected by the instinctive reaction to turn away and blink.
2. Observe the regulations valid in the respective country where the experiment is being performed, e.g. Germany's health and safety regulations BGV B2 "Laser radiation", and any stipulations set by the relevant ministry.
3. **Do not stare into the laser beam.** Accidental short exposures are normally limited by the natural blink reflex, but intentional staring into the beam can cause eye injury.
4. **Do not aim the laser beam at people or animals.** Always keep the beam path below eye level and safely terminate it with a beam stop.
5. **Avoid reflective surfaces.** Reflections may be as hazardous as direct exposure.
6. **Do not use optical instruments (e.g. magnifier, binoculars, lens) to view the laser beam.**

Use in Schools

7. Operation is permitted **only under the supervision of a qualified teacher or trained supervisor.**
8. **Before experiments begin**, students must be informed about potential eye hazards and correct handling.
9. The laser must be **stored in a locked place when not in use** and must not be operated by students without supervision.
10. The device must **not be used for play, competitions or games.**

Handling and Operation

11. Safe operation of the laser is guaranteed, provided it is used correctly. However, there is no guarantee of safety if the equipment is used in an inappropriate or careless manner. If it is deemed that the equipment can no longer be operated without risk (e.g. visible damage has occurred), the laser should be switched off immediately and secured against any unintended use.
12. Before putting the equipment into operation, check for any signs of damage. In the event of any malfunction or visible damage, turn off the laser and put it away so that it cannot be used unintentionally.
13. Any modification, opening or alteration of the device is prohibited.
14. Always keep the device switched off during transport and, if available, use the **protective cap over the output aperture.**
15. Do not conduct experiments where the beam runs freely across long distances. Always set up experiments so that the beam is safely terminated by a **beam stop / screen.**
16. For long-term setups or continuous use, **mark or block off the laser area** so the beam path is clearly visible.

Maintenance and Storage

17. Clean only with a soft, dry cloth; do not use solvents or chemicals.
18. Protect the device from dust, moisture and high temperatures.
19. Use the specified power supply only. Do not use damaged cables or plugs.

Emergency Procedure

20. If eye exposure is suspected, immediately inform the supervising teacher and seek ophthalmological examination.
21. Any irregularities (lose parts, damaged lens, defective cables) must be reported at once; the device must not be used until inspected.

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Specifications

EN

Laser Safety	Class II
Power	< 1 mW
Wavelength	532 nm ± 10 nm
Operating Temperature	-30°C to +50°C
Laser Type	Continuous Wave Laser
Laser Technology	DPSS laser (diode-pumped solid-state laser)
Power Supply	Primary: 100V – 240V, 50/60Hz Secondary: 5 V DC, 1000mA
Environmental Conditions	Indoor use only IP20, 5 – 35°C

Notizen / Notes