


INTRAMUSCULAR INJECTION SIMULATOR – UPPER ARM

This high-quality I.M. injection skill trainer represents a right upper arm with all important anatomical **palpable landmarks** such as **acromion and humerus**.

The realistic anatomy allows for placing **correct intramuscular injections** in an exceptionally graphic way.



➤ INTRAMUSCULAR INJECTION CAN BE PRACTICED REALISTICALLY WITHOUT A REAL PATIENT!



Realistic anatomical landmarks such as acromion and humerus



Providing both training and control modes for an immediate or subsequent performance check



Audio-visual feedback on the injection performance regarding site, depth and bone contact



Exchangeable lifelike silicone skin

Delivery content:

- P55 Injection Simulator (upper arm) with humerus stump, (foam) muscle and silicone skin preassembled
- Replacement silicone skin
- Disposable 5 ml syringe
- Injection cannula for intramuscular injections
- AA batteries

Dimensions and weight: 32 x 14 x 12 cm, 0,9 kg

1009840 Intramuscular Injection Simulator – Upper Arm



01/2021

 MORE INFORMATION AVAILABLE ONLINE AT 3BSCIENTIFIC.COM

