



BOVINE VENIPUNCTURE & EPIDURAL TAIL

The VSI Bovine Venipuncture & Epidural Tail is designed to seamlessly integrate with existing VSI Bovine Dystocia and Theriogenology simulators or function independently on a tabletop stand. Made from durable, lifelike silicone, it provides a realistic texture and flexibility that enhances hands-on learning. It comes with two tail options, modeled after Holstein and Hereford breeds. This training tool is ideal for practicing key veterinary procedures such as venipuncture and epidural placement, allowing students to develop precision and confidence in a controlled, realistic environment.

TECHNICAL SPECIFICATIONS:

- Silicone tail
- Compatible with VSI's Bovine Dystocia and Bovine Theriogenology models
- Includes tabletop stand for independent use
- Full articulation
- Venipuncture and Epidural Functions
- Vacuum draw when epidural is performed
- Includes VSI synthetic blood
- Durable Hide covering
- Features palpable replica vertebrae



DAY ONE SKILLS & PRODUCT BENEFITS:

- Venipuncture Training: Practice accurate blood draws using the realistic silicone tail.
- Epidural Placement: Train in administering epidurals in a controlled, life like setup.



Hereford Epidural/Venipuncture Tail – Black
1025526

Hereford Epidural/Venipuncture Tail – Brown
1025527

Holstein Epidural/Venipuncture Tail
1025536

REPLACEABLE PARTS:

- | | |
|--|---------|
| • Bovine Venipuncture Tail Veins (Set of 10) | 1026408 |
| • Bovine Tail Stand | 1025538 |
| • Neck Venipuncture Fluid | 1025547 |
| • Bovine Epidural Tubing (Set of 10) | 1026409 |
| • Slow Release Valve | 1026412 |

06/2025 9001463

IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



3B Scientific Group

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · [3bscientific.com](https://www.3bscientific.com)