




CANINE BASIC LIFE SUPPORT SIMULATOR

The VSI Canine Basic Life Support simulator features an articulated skeleton with a ribcage, specifically designed for effective chest compression practice. It includes a feedback monitor that provides real-time information on compression speed and depth, helping students refine their techniques efficiently.

The model also offers realistic representations of the trachea and esophagus, enabling practice in intubation and nasogastric tube placement. Additionally, it includes jugular, cephalic, and saphenous veins for practicing venipuncture and catheterization, making it a comprehensive resource for developing critical clinical skills. Ideal for practicing restraint, positioning, and bandaging techniques, this model also features representations of the ear and simulated debris, allowing students to learn proper ear cleaning methods.



TECHNICAL SPECIFICATIONS

- Canine BLS Simulator
- 2x AA Batteries
- Ear Wax
- Cradle
- Blood Bag
- Veins
- Carry Bag

DAY ONE SKILLS & PRODUCT BENEFITS

- Fully articulated steel skeleton with ribcage, and separated toes
- Representation of trachea and esophagus which allows for practical intubation training
- Representation of nasal passage to allow for nasogastric intubation training
- Representation of jugular, cephalic, and saphenous veins to practice venipuncture and catheterization in front and hint limbs
- Chest compression function with feedback monitor which provides feedback on speed and depth of compressions
- Representation of ear and simulated debris to practice ear cleaning techniques
- Ideal for practicing restraint and bandaging techniques

REPLACEABLE PARTS

- Venipuncture Tubing 1025634
- Hide Covering 1025624
- Venipuncture Fluid Bag 1025663
- Tongue
- Laryngopharynx Component
- Ear Capsule (Left & Right)
- Compression Monitor

IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



3B Scientific Group

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · 3bscientific.com