



The frame has been streamlined to minimize the wearer's profile, increasing the Hyper-Realism of any training scenario. The TCCC Cut Suit is fully user-customizable and user-repairable, reducing costly repairs and long-term recurring expenses.



UNIQUE FEATURES:

- The system can be worn during intense physical scenarios at the Point of Injury (POI)
- The system weighs approximately 20 lbs (9,07 kg)
- Clothing, uniform, body armor, and equipment are usable
- The system allows for interaction with a live human during the emergency assessment and treatment process
- The skin is user repairable
- User created and customizable wounds

PROCEDURES CURRENTLY AVAILABLE ON THE TCCC CUT SUIT:

- Extremity hemorrhage control with tourniquet application, arterial ligation / clamping, internal compression (wound void packing)
- Surgical Cricothyroidotomy
- Bilateral anterior and axillary chest needle thoracocentesis
- Bilateral surgical chest-tube thoracotomy
- Suturing and stapling of skin
- Peripheral IV access
- Sternal intraosseous IV access
- Peripheral IV medication administration



Field expedient cricothyroidotomy



Chest needle thoracocentesis



Occlusive dressing



Chest tube



Hemorrhage control - arm



Hemorrhage control - leg

TCCC/EMS Cut Suit - Light Skin | 1026000
TCCC/EMS Cut Suit - Dark Skin | 1025995

IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



3B Scientific Group

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · 3bscientific.com