


NEW



AIRWAY MANAGEMENT TRAINER P73

EXTENTION TO ATLAS – NEW 3B SCIENTIFIC AIRWAY TRAINER

Introducing an anatomically accurate, realistic model designed for comprehensive training in airway management techniques. The model replicates the weight and feel of a human head for a realistic training experience. It facilitates health-care professionals and students to practice essential skills like endotracheal intubation and the use of supraglottic airway devices, ensuring proficiency in real-world situations.

ACHIEVE PROFICIENCY IN AIRWAY MANAGEMENT TECHNIQUES

Ideal for healthcare educators, medical students, and simulation centers, this trainer offers an essential solution for mastering airway management in a controlled, safe and effective learning environment.



Nasotracheal



LT Insertion



BVM Ventilation



i-gel Insertion

REALISTIC AIRWAYS

- ✓ Tongue
- ✓ Nasopharynx
- ✓ Oropharynx
- ✓ Epiglottis
- ✓ Vallecula
- ✓ Vocal cords
- ✓ Thyroid cartilage
- ✓ Esophagus
- ✓ Trachea
- ✓ Lung
- ✓ Stomach



MORE FEATURES

- Oral and nasal intubation using standard basic and advanced airway management devices:
 - ✓ OPA, ✓ NPA, ✓ i-gel, ✓ LMA, ✓ LT, ✓ ETT
- Right Main stem Intubation
- Esophageal Intubation
- Gastric Insufflation
- Teeth Damage Indication
The upper tooth design reveals potential tooth damage risks during laryngoscope use, helping users refine their technique and gain patient safety.
- Three different wide pupils (Normal/Constricted/Dilated)



CONFIGURATIONS

- Skin color:
 - Light skin 1026007
 - Dark skin 1026076
- Optional 4th Order Segmental Bronchi 1026077
- Delivery Contents: Airway Management Trainer, Lubricant, Carry Bag



IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE



3B Scientific Group

Headquarters · Ludwig-Erhard-Strasse 20 · 20459 Hamburg · Germany · [3bscientific.com](https://www.3bscientific.com)