



Advance Bovine Theriogenology Model

Advanced Training for Artificial Insemination & Reproductive Procedures

The Hereford Advanced Bovine Theriogenology Model is a premium, life-like training system designed to provide students with realistic, hands-on experience in bovine reproductive health, palpation, and artificial insemination. Built for durability and anatomical accuracy, this model prepares learners for real clinical and field scenarios in veterinary and animal science careers.

PRODUCT FEATURES:

- Rear section constructed from durable epoxy/fiberglass
- Life-like pelvis and soft perineum palpation panel
- Inflatable vinyl rectum for realistic palpation resistance
- Flexible tail for procedural positioning
- Integrated venipuncture and epidural tail
- Functional milkable udders with integrated milk tank supporting California Milk Test procedures
- Theriogenology Uterus Set with 5 interchangeable models, featuring detailed cervix, ovaries, and broad ligament anatomy with realistic ovarian structures such as follicles and corpus luteum:
 - Open Cow – Ideal for artificial insemination training
 - Open Heifer – Ideal for artificial insemination training
 - 45 Days Pregnant – Enlarged horn with amniotic vesicle
 - 60 Days Pregnant – Uterine slip and detectable fetus
 - 90 Days Pregnant – Uterine slip, fetus, and placentomes

Hereford Advanced Bovine Theriogenology Model,
Black | 1026384
Hereford Advanced Bovine Theriogenology Model,
Brown | 1026388
Hereford Advanced Bovine Theriogenology Model,
Brown White | 1026389



DAY ONE SKILLS & PRODUCT BENEFITS

- Rectal palpation and reproductive anatomy identification
- Artificial insemination technique for cows and heifers
- Pregnancy diagnosis at multiple gestational stages
- Ovarian evaluation, including follicular and luteal structures
- Dairy testing procedures using milkable udders
- Large-animal venipuncture and epidural administration practice
- Anatomical orientation of the cervix, uterus, ovaries, and broad ligament
- Durable construction designed for repeated instructional use

REPLACEABLE PARTS:

Bovine Ovary Set	1025523
Bovine Theriogenology Uterus Set (Set of 5)	1025499
Hereford Palpation Panel Assembly Black	1025492
Hereford Palpation Panel Assembly Brown	1025493
Hereford Teats Black	1026404
Hereford Teats	1025495
Hereford Milk Tank	1026474
Bovine Epidural Tubing (Set of 10)	1026409
Bovine Venipuncture Tail Veins (Set of 10)	1026408
Neck Venipuncture Fluid	1025547



02/2026

IF YOU HAVE ANY QUESTIONS OR NEED SUPPORT, PLEASE SCAN HERE

