

ON-THE-GO RECORDING AND DEBRIEFING SYSTEM

NORDIC
CoReCMOBILE



#SmartAndEasyLearning

PORTABLE & WIRELESS RECORDING SYSTEM

Introducing the Nordic CoRecMobile – a cutting-edge, portable recording system designed for capturing a wide-range of learning experiences with ease, anywhere, anytime.


Packaged in a convenient wheeled protective case, it's perfect for any environment. As part of the Nordic Simulators product family, the CoRecMobile is powered by the intuitive Nordic CoRecorder software, making it user-friendly and efficient. The system comes equipped with three high-definition wireless cameras, all featuring magnetic and clamp fastenings for quick setup.

Developed in collaboration with simulation experts, Nordic CoRecMobile stands out for its high quality, versatility, and ease of use. Turn any space into a simulation studio in ten minutes.

Key Features include Wireless Pan-Tilt-Zoom cameras with built-in microphones and software controlled AI human tracking. The cameras have internal batteries for complete mobility - take your simulation exercises out of the training room and On-the-Go.

Elevate your simulation training with Nordic CoRecMobile – designed for professionals, by professionals.

WHAT'S INCLUDED



PLUG AND PLAY

- 1 Extremely easy to use intuitive browser-based user interface
- 2 Cover multiple rooms at once with three synchronized cameras
- 3 Control your camera views and AI tracking directly from the software
- 4 View, record and debrief multiple cameras and audio feeds simultaneously
- 5 Insert customizable chapter marks and comments into the recordings
- 6 For easy debriefing - edit and modify your recordings
- 7 Export your recordings for storage and sharing
- 8 Capture and record any HDMI source as one of the cameras



CONTACT US

office@nordicsimulators.fi

+358 (03) 589 8588

Vuoripojankatu 16
15210 Lahti Finland



www.nordicsimulators.com



facebook.com/NordicSimulators



instagram.com/nordicsimulators



x.com/nordicsimulator